

Create and use custom devices in ZEN Blue using the Arduino microcontroller

- How to leverage the full potential of ZEN Blue
and 3rd party image processing (or hardware)
in your image processing and analysis workflows
(and image acquisition)


Usual choices for development of home-made custom devices

- Write a new software from scratch
- Labview
- Matlab (ScanImage)
- Micro-Manager

- Proprietary software (eg ZEN)

Lumalum, a root imaging robot

with H    ne Javot and Michel Philibert (CEA)



java com layer
mm bsh
mm plugin
ij macros
...


Simplifying MM device interaction

```
run("MM MacroExtensions");
```

```
for ( i = 0 ; i < 50 ; i = i+2 ) {  
    Ext.moveRelativeXYZ( 0 , 0 , i );  
    Ext.snap( );  
}
```

```
run ("Images to Stack");  
run ("EDF Easy ", "quality='0' topology='1'");
```

MM design model




<https://micro-manager.org/wiki/Micro-Manager%20Project%20Overview>



Zeiss MTB design model

1.1 Components of the MTB2011

Carl Zeiss MicroToolbox 2011




Elements of device interaction



Target: Axiozoom under Zen Blue
Ability to control external light from Zen




The Zeiss *.mtb device model & Micro Tool Box Configuration




```
<MTB2004Export Version="2.7.0.5">
  <MTBConfiguration Name="New Configuration1" MyChildId="0">
    <Devices>
      <Device Class="MTBCtrlDevice,MTBKernel" Name="Wega.Stand" External="false"
        PortType="None" PortNo="1" PortBaudAddr="57600" RTSCTS="false" CanNode="25"
        MyChildId="1">...</Device>
      <Device Class="MTBCtrlElectronicModuleDeviceREO,MTBKernelREO"
        Name="Wega_Gateway_V16.EMS3" External="false" PortType="USB" PortNo="4111"
        PortBaudAddr="" RTSCTS="false" CANConnection="" CanNode="40" MyChildId="11600">...
        </Device>
      <Device Class="MTBCtrlDeviceSycop3,MTBKernelWEGA" Name="Wega.Main" External="false"
        PortType="USB" PortNo="4111" PortBaudAddr="" RTSCTS="false" CanNode="25" MyChildId="0">
        <SimulationMode>true</SimulationMode>
        <SimulationClass>ZEISS.MTB.Simulation.MTBSimDeviceWegaMain</SimulationClass>
        <TechReport>ReportWegaMain</TechReport>
        <Components>
          <Component Class="MTBCtrlMicroscopeManagerSycop3,MTBKernelREO"
            Name="REO.MicroscopeManager" MTBId="MTBMicroscopeManager" CANId="00"
            MyChildId="11024">
              <LightManagerModes>Off</LightManagerModes>
            </Component>
        </Components>
      </Device>
      <Device Class="MTBCtrlOpticsDeviceWEGA,MTBKernelWEGA" Name="Wega.OpticsDevice_mot"
        External="false" PortType="USB" PortNo="4111" PortBaudAddr="" RTSCTS="false"
        CanNode="32" MyChildId="0">...</Device>
      <Device Class="MTBCtrlFocusDeviceREO,MTBKernelREO" Name="Wega.FocusDevice_mot"
        External="false" PortType="USB" PortNo="4111" PortBaudAddr="" RTSCTS="false"
        CanNode="36" CANConnection="" MyChildId="11500">...</Device>
      <Device Class="MTBCtrlOpticsDeviceREO,MTBKernelREO" Name="Wega.FLTube_mot"
        External="false" PortType="USB" PortNo="4111" PortBaudAddr="" RTSCTS="false"
        CanNode="102" CANConnection="" MyChildId="11502">...</Device>
      <Device Active="false" Name="ApoTome.none" MyChildId="11604"/>
      <Device Class="MTBCtrlCan29HXF120Device,MTBKernelCan29Changers" Name="HXP120.Device"
        External="false" PortType="USB" PortNo="4099" PortBaudAddr="" RTSCTS="false"
        CANConnection="" CanNode="54" MyChildId="11601">...</Device>
      <Device Active="false" Name="Stage.none" MyChildId="11602"/>
      <Device Active="false" Name="Incubation.AtmoTemp_none" MyChildId="14000"/>
      <Device Active="false" Name="Incubation.FlowTemp_none" MyChildId="14001"/>
      <Device Active="false" Name="Linkam.none" DisplayName="No Linkam Heating Stage"
        DisplayName-de="Kein Linkam Heiztisch" MyChildId="9914"/>
      <Device Active="false" Name="IOCard.NoIOCard" DisplayName="No IOCard" DisplayName-
        de="Keine IO-Karte" External="true" PortType="RS232" PortNo="2" PortBaudAddr="9600"
        RTSCTS="true" MyChildId="25100"/>
      <Device Active="false" Name="NoFilterWheelsShutters" DisplayName="No Filter Wheels and
        Shutters" DisplayName-de="Keine Filterräder und Verschlüsse" External="true"
        PortType="RS232" PortNo="2" PortBaudAddr="9600" RTSCTS="true" MyChildId="30050"/>
    </Devices>
  </MTBConfiguration>
</MTB2004Export>
```

Chosen device : Sutter Lambda 10-3



Lambda 10-3
Optical Filter Changer

Emulating the Sutter protocol




COM2 (Arduino Uno) - Serial Port Monitor - [Table view]

| # | Time | Function | Dire... | Status | Data | Data (chars) |
|----|---------------------|---|---------|----------------|--------------------|---------------|
| 5 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_LINE_CONTROL) | UP | STATUS_SUCCESS | 00 00 00 | ... |
| 6 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_CHARS) | DOWN | | | |
| 7 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_CHARS) | UP | STATUS_SUCCESS | 00 00 00 00 00 ... | |
| 8 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_HANDFLOW) | DOWN | | | |
| 9 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_HANDFLOW) | UP | STATUS_SUCCESS | 00 00 00 00 40 ... |@..... |
| 10 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_BAUD_RATE) | DOWN | | 80 25 00 00 | €%.. |
| 11 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_BAUD_RATE) | UP | STATUS_SUCCESS | | |
| 12 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_RTS) | DOWN | | | |
| 13 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_RTS) | UP | STATUS_SUCCESS | | |
| 14 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_DTR) | DOWN | | | |
| 15 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_DTR) | UP | STATUS_SUCCESS | | |
| 16 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_LINE_CONTROL) | DOWN | | 00 00 08 | ... |
| 17 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_LINE_CONTROL) | UP | STATUS_SUCCESS | | |
| 18 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_CHARS) | DOWN | | 00 00 00 00 11 ... | |
| 19 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_CHARS) | UP | STATUS_SUCCESS | | |
| 20 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_HANDFLOW) | DOWN | | 01 00 00 00 40 ... |@...d..d. |
| 21 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_HANDFLOW) | UP | STATUS_SUCCESS | | |
| 22 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_TIMEOUTS) | DOWN | | e8 03 00 00 00 ... | è.....è..... |
| 23 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_TIMEOUTS) | UP | STATUS_SUCCESS | | |
| 24 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_QUEUE_SIZE) | DOWN | | e8 03 00 00 e8 ... | è...è... |
| 25 | 26/07/2019 12:15:21 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_QUEUE_SIZE) | UP | STATUS_SUCCESS | | |
| 26 | 26/07/2019 12:15:21 | IRP_MJ_READ | DOWN | | | |
| 27 | 26/07/2019 12:15:21 | IRP_MJ_WRITE | DOWN | | ee | ì |
| 28 | 26/07/2019 12:15:21 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | ee | ì |
| 29 | 26/07/2019 12:15:21 | IRP_MJ_READ | UP | STATUS_SUCCESS | ee | ì |
| 30 | 26/07/2019 12:15:21 | IRP_MJ_READ | DOWN | | | |
| 31 | 26/07/2019 12:15:21 | IRP_MJ_READ | UP | STATUS_SUCCESS | 0d | . |
| 32 | 26/07/2019 12:15:21 | IRP_MJ_READ | DOWN | | | |
| 33 | 26/07/2019 12:15:22 | IRP_MJ_READ | UP | STATUS_TIMEOUT | | |
| 34 | 26/07/2019 12:15:22 | IRP_MJ_READ | DOWN | | | |
| 35 | 26/07/2019 12:15:22 | IRP_MJ_WRITE | DOWN | | fd | ý |
| 36 | 26/07/2019 12:15:22 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | fd | ý |
| 37 | 26/07/2019 12:15:22 | IRP_MJ_READ | UP | STATUS_SUCCESS | fd | ý |
| 38 | 26/07/2019 12:15:22 | IRP_MJ_READ | DOWN | | | |
| 39 | 26/07/2019 12:15:22 | IRP_MJ_READ | UP | STATUS_SUCCESS | 31 | 1 |


For help press F1

IRP: 51 Read: 33 Written: 2



ZEN <-> RS232 <-> LB10-3

https://www.sutter.com/manuals/LB10-3_OpMan.pdf
<https://www.arduino.cc/>
<https://www.eltima.com/products/serial-port-monitor/>




J.Mutterer, MIAP OAD workshop, Freiburg

MTB connection test

Filter Wheels and Shutters(Lambda 10-3)

Connection

Type: Port:



| | | | | | | | | |
|-----|---------------------|---|------|----------------|--------------------|--------------------|----|------|
| 102 | 26/07/2019 12:23:46 | IRP_MJ_CREATE | DOWN | | | C:\Program File... | | COM2 |
| 103 | 26/07/2019 12:23:46 | IRP_MJ_CREATE | UP | STATUS_SUCCESS | | C:\Program File... | | COM2 |
| 104 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_BAUD_RATE) | DOWN | | | | | COM2 |
| 105 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_BAUD_RATE) | UP | STATUS_SUCCESS | 80 25 00 00 | €%.. | 4 | COM2 |
| 106 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_LINE_CONTROL) | DOWN | | | | | COM2 |
| 107 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_LINE_CONTROL) | UP | STATUS_SUCCESS | 00 00 08 | ... | 3 | COM2 |
| 108 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_CHARS) | DOWN | | | | | COM2 |
| 109 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_CHARS) | UP | STATUS_SUCCESS | 00 00 00 00 11 ... | | 6 | COM2 |
| 110 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_HANDFLOW) | DOWN | | | | | COM2 |
| 111 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_GET_HANDFLOW) | UP | STATUS_SUCCESS | 01 00 00 00 40 ... |@...d...d... | 16 | COM2 |
| 112 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_BAUD_RATE) | DOWN | | 80 25 00 00 | €%.. | 4 | COM2 |
| 113 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_BAUD_RATE) | UP | STATUS_SUCCESS | | | | COM2 |
| 114 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_RTS) | DOWN | | | | | COM2 |
| 115 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_RTS) | UP | STATUS_SUCCESS | | | | COM2 |
| 116 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_DTR) | DOWN | | | | | COM2 |
| 117 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_DTR) | UP | STATUS_SUCCESS | | | | COM2 |
| 118 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_LINE_CONTROL) | DOWN | | 00 00 08 | ... | 3 | COM2 |
| 119 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_LINE_CONTROL) | UP | STATUS_SUCCESS | | | | COM2 |
| 120 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_CHARS) | DOWN | | 00 00 00 00 11 ... | | 6 | COM2 |
| 121 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_CHARS) | UP | STATUS_SUCCESS | | | | COM2 |
| 122 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_HANDFLOW) | DOWN | | 01 00 00 00 40 ... |@...d...d... | 16 | COM2 |
| 123 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_HANDFLOW) | UP | STATUS_SUCCESS | | | | COM2 |
| 124 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_TIMEOUTS) | DOWN | | e8 03 00 00 00 ... | è.....è..... | 20 | COM2 |
| 125 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_TIMEOUTS) | UP | STATUS_SUCCESS | | | | COM2 |
| 126 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_QUEUE_SIZE) | DOWN | | e8 03 00 00 e8 ... | è...è... | 8 | COM2 |
| 127 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_QUEUE_SIZE) | UP | STATUS_SUCCESS | | | | COM2 |
| 128 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | 1 | COM2 |
| 129 | 26/07/2019 12:23:46 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | ee | ì | 1 | COM2 |
| 130 | 26/07/2019 12:23:46 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | ee | ì | 1 | COM2 |

Protocol standard : not standard

| | | |
|------------------------|-----------------------|-------------------------|
| Transfer to On Line | 238 EE 11101110 | Puts controller on-line |
|------------------------|-----------------------|-------------------------|

5.6.2 Command Transmission Protocol


The Lambda 10-3 does not use any of the standard protocols commonly used for serial line or USB communications between computers or between a computer and a peripheral device. The Lambda 10-3 controller, however, does generate a primitive form of protocol of which control software running on the remote computer can, and should, make use. This protocol consists of two main components: “confirmation command echo” and “command completion indicator”. This command transmission protocol is used in the same manner for RS-232 serial port, USB, and parallel port connections.

5.6.2.1 Confirmation Command Echo

When the host computer sends a command to the Lambda 10-3, each byte received is immediately echoed back to the host computer. This echoing back of each sent byte is a confirmation that the byte has been received and will be acted upon shortly. A short period after the last byte of a command has been echoed back, the Lambda 10-3 sends a confirmation byte (described next) that indicates the operation associated with the command completed.

5.6.2.2 Command Completion Indicator

When the Lambda 10-3 completes the operation associated with the command it has just received, it transmits back to the host computer a byte value of 13 decimal (0D hexadecimal, 00001101 binary). This byte value corresponds to an ASCII carriage return (often abbreviated as “CR”).



```
byte cmd ;

void setup() {
  Serial.begin(9600);
  pinMode(13,OUTPUT);
  digitalWrite(13,LOW);
}

void loop()
{
  if (Serial.available()) {
    cmd = Serial.read();
    processCommand(cmd);
  }
}

void processCommand(byte s) {
  if (s==0xee) {
    Serial.write(s); // echo ONLINE cmd
    Serial.write(0x0d);
  }
}
```

Configuration command **0xFD** : simple echo not enough

Filter Wheels and Shutters(Lambda 10-3)

Connection

Type:

RS232

Port:



2

Check Connection

✖

| | | | | | | | | | |
|-----|---------------------|---|------|----------------|----|---|---|---|------|
| 127 | 26/07/2019 12:23:46 | IRP_MJ_DEVICE_CONTROL (IOCTL_SERIAL_SET_QUEUE_SIZE) | UP | STATUS_SUCCESS | | | | | COM2 |
| 128 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 129 | 26/07/2019 12:23:46 | IRP_MJ_WRITE | DOWN | | ee | i | 1 | 1 | COM2 |
| 130 | 26/07/2019 12:23:46 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | ee | i | 1 | | COM2 |
| 131 | 26/07/2019 12:23:46 | IRP_MJ_READ | UP | STATUS_SUCCESS | ee | i | 1 | | COM2 |
| 132 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 133 | 26/07/2019 12:23:46 | IRP_MJ_READ | UP | STATUS_SUCCESS | 0d | . | 1 | | COM2 |
| 134 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 135 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_TIMEOUT | | | | | COM2 |
| 136 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 137 | 26/07/2019 12:23:47 | IRP_MJ_WRITE | DOWN | | fd | y | 1 | 1 | COM2 |
| 138 | 26/07/2019 12:23:47 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | fd | y | 1 | | COM2 |
| 139 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | fd | y | 1 | | COM2 |
| 140 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 141 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | 0d | . | 1 | | COM2 |
| 142 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 143 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_TIMEOUT | | | | | COM2 |
| 144 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |

Just the model **(10-3)** is enough for MTB test check, but not for ZEN ! **x**



| | | | | | | | | | |
|-----|---------------------|--------------|------|----------------|----|---|---|---|------|
| 128 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 129 | 26/07/2019 12:23:46 | IRP_MJ_WRITE | DOWN | | ee | ì | | 1 | COM2 |
| 130 | 26/07/2019 12:23:46 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | ee | ì | | 1 | COM2 |
| 131 | 26/07/2019 12:23:46 | IRP_MJ_READ | UP | STATUS_SUCCESS | ee | ì | | 1 | COM2 |
| 132 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 133 | 26/07/2019 12:23:46 | IRP_MJ_READ | UP | STATUS_SUCCESS | 0d | . | | 1 | COM2 |
| 134 | 26/07/2019 12:23:46 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 135 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_TIMEOUT | | | | | COM2 |
| 136 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 137 | 26/07/2019 12:23:47 | IRP_MJ_WRITE | DOWN | | fd | ý | | 1 | COM2 |
| 138 | 26/07/2019 12:23:47 | IRP_MJ_WRITE | UP | STATUS_SUCCESS | fd | ý | | 1 | COM2 |
| 139 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | fd | ý | | 1 | COM2 |
| 140 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |
| 141 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | 31 | . | 1 | 1 | COM2 |
| 142 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | . | | 1 | COM2 |
| 143 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | 30 | . | 0 | 1 | COM2 |
| 144 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | . | | 1 | COM2 |
| 145 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | 2d | . | - | 1 | COM2 |
| 146 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | . | | 1 | COM2 |
| 147 | 26/07/2019 12:23:47 | IRP_MJ_READ | UP | STATUS_SUCCESS | 33 | . | 3 | 1 | COM2 |
| 148 | 26/07/2019 12:23:47 | IRP_MJ_READ | DOWN | | | | | 1 | COM2 |

Fully functional reply to configuration command 0xFD


5.4.9 Get Controller Type and Configuration

This command is used to obtain information about the controller as to its model and configuration. The following table shows the type of information returned when issuing this command

Table 5-8. "Get Controller Type and Configuration" command return codes and data.

| Total Num. Bytes | Description | | | |
|------------------|---------------------|------------|-----------------|-----------------|
| | Category | Num. Bytes | Possible Values | |
| | | | ASCII String | Meaning |
| 29 | Controller Type | 4 | 10-3 | Lambda 10-3 |
| | Filter Wheel A Type | 5 | WA-25 | 25mm |
| | | | WA-32 | 32mm |
| | | | WA-HS | High Speed |
| | | | WA-BD | Belt Driver |
| | | | WA-NC | Not Connected |
| | | | WA-ER | Error |
| | Filter Wheel B Type | 5 | WB-25 | 25mm |
| | | | WB-32 | 32mm |
| | | | WB-HS | High Speed |
| | | | WB-BD | Belt Driver |
| | | | WB-NC | Not Connected |
| | | | WB-ER | Error |
| | Filter Wheel C Type | 5 | WB-25 | 25mm |
| | | | WB-32 | 32mm |
| | | | WB-HS | High Speed |
| | | | WB-BD | Belt Driver |
| | | | WB-NC | Not Connected |
| | | | WB-ER | Error |
| | Shutter A Type | 5 | SA-IQ | SmartShutter |
| | | | SA-VS | Vincent Shutter |
| | Shutter B Type | 5 | SB-IQ | SmartShutter |
| | | | SA-VS | Vincent Shutter |

```
void processCommand(byte s) {
  if (s == 0xee) {
    Serial.write(s); // echo ONLINE cmd
    Serial.write(0x0d);
  } else if (s == 0xfd) {
    Serial.write(s); // echo CONFIG cmd
    Serial.print("10-3WA-BDWB-NCWC-NCSA-VSSB-VS");
    Serial.write(0x0d);
  }
}
```



Full reply to the status command 0xCC

5.4.2 Status

The Status command is used to return information about the filter wheel state (if installed) and the state of one or two shutters (if, and as, installed). The following table describes the type of data and information returned by the Status command.

Table 5-7. Status command return codes and data.




| Order | Num. of bytes | Category | Sub Category | Value (Decimal, hexadecimal, & binary) | Description |
|-------|---------------|------------------------------|---------------------------|---|---|
| 1 | 1 | Command echo | | 204 CC 11001100 | The Status command byte code echoed back. |
| 2 | 1 | Filter Wheel A status | Wheel, Speed and Position | Encoded command byte: 0 – 121 00 – 79 00000000 – 01111001 * | Filter Wheel (Bit 7): 0 |
| | | | | | Speed (Bits 6, 5, & 4): 0 – 7. |
| | | | | | Position (Bits 3, 2, 1, & 0): 0 – 9 |
| 5 | 1 | Shutter A open/-closed state | Open | 170 AA 10101010 | Shutter A is in the open state. |
| | | | Open state is conditional | 171 AB 10101011 | The open state of Shutter A is conditional upon the movement of the Filter Wheel A. |
| | | | Closed | 172 AC 10101100 | Shutter A is in the closed state. |

```
void processCommand(byte s) {
  if (s == 0xee) {
    Serial.write(s); // echo ONLINE cmd
    Serial.write(0x0d);
  } else if (s == 0xfd) {
    Serial.write(s); // echo CONFIG cmd
    Serial.print("10-3WA-BDWB-NCWC-NCSA-VSSB-VS");
    Serial.write(0x0d);
    Serial.flush();
  } else if (s == 0xcc) {
    Serial.write(s); // echo STATUS cmd
    Serial.write(0x00); // FWA
    Serial.write(0x80); // FWB
    Serial.write(0xfc); // FWC
    Serial.write(0x00);
    Serial.write(0xac); // SA closed
    Serial.write(0xbc); // SB closed
    Serial.write(0xdc); // SA SmartShutter Fast mode
    Serial.write(0x01);
    Serial.write(0xdc); // SB SmartShutter Fast mode
    Serial.write(0x02);
    Serial.write(0x0d); // confirmation
    Serial.flush();
  }
}
```




Functional Shutter

| | | | | | | |
|------|----------------|----|---|---|---|------|
| DOWN | | aa | a | 1 | 1 | COM2 |
| UP | STATUS_SUCCESS | aa | a | 1 | | COM2 |
| UP | STATUS_SUCCESS | aa | a | 1 | | COM2 |
| DOWN | | | | | 1 | COM2 |
| UP | STATUS_SUCCESS | 0d | . | 1 | | COM2 |
| DOWN | | | | | 1 | COM2 |
| UP | STATUS_TIMEOUT | | | | | COM2 |
| DOWN | | | | | 1 | COM2 |
| DOWN | | ac | ┐ | 1 | 1 | COM2 |
| UP | STATUS_SUCCESS | ac | ┐ | 1 | | COM2 |
| UP | STATUS_SUCCESS | ac | ┐ | 1 | | COM2 |
| DOWN | | | | | 1 | COM2 |
| UP | STATUS_SUCCESS | 0d | . | 1 | | COM2 |

```
} else if (s==0xaa) {  
  Serial.write(s);          // echo SA OPEN  
  Serial.write(0x0d);  
  digitalWrite(13,HIGH);  
} else if (s==0xac) {  
  Serial.write(s);          // echo SA CLOSE  
  Serial.write(0x0d);  
  digitalWrite(13,LOW);  
} else {  
  Serial.write(s); // echo any other case  
  Serial.write(0x0d);  
}
```



Construction of the objective support



Levels of Zen Blue 2.3 integration



Microscope Control

50x

100% Came... Test Camera

Test Camer... 100% Obse...

none ApoTome Closed Closed 0 %

5x / 0.16

Bright Field 0 % 0 V

Microscope... TL Vert Closed ApoTome P...

Favorites

Configure...

on off

Experiment Feedback

Show All

Edit Feedback Script...

Select script runtime conditions

Free Run Synchronized

Script runs freely

Acquisition Analysis HD Writing

Script Run*

* when triggered by used observable inside the LoopScript

Macro Editor

File Edit Record Debug Help

New Macro Record Run Debug Step Over

demo toulouse Macro-02

```
1 import time
2
3 def leds(name,t):
4     ... led.Load(name)
5     ... Zen.Devices.ApplyHardwareSetting(led)
6     ... time.sleep(t)
7
8 led = ZenHardwareSetting()
9 leds("on",0)
10 leds("blue",1)
11 leds("white",2)
12 leds("red",1)
13 leds("off",0)
```

User Documents

demo toulouse

Macro-01

Macro-02

Link to this presentation:

<https://docs.google.com/presentation/d/1BIDqrpJWmzF-OCsNnKX5SfcnYt5UuXP2ksjkCaLP0ZQ/edit?usp=sharing>

Link to the github repository with example code for emulated devices:

https://github.com/mutterer/MM_CustomArduino